



# Skedco® Kiddie Litter Insert with Haz-Mat/Hospital Sked®

**Weight:** 3 pounds (Insert Alone), 15 pounds (Insert with HMH Sked)

**Dimension:** 36 inch wide x 96 inch long

**Materials:** Insert: 420 Denier packcloth nylon fabric.  
Pockets: 13.5 oz. Polyester double sided arctic fleece  
HMH Sked: Medium density E-Z glide non-absorbing polyethylene plastic. 0.105 inch thick.

**Product Number:** Kiddie Litter Insert: SK-255  
Kiddie Litter Insert with HMH Sked: SK-2550



The Kiddie Litter was developed for emergency evacuation of babies from the nursery or other locations in hospitals in a safer manner than anything available anywhere else. It is an insert that fits in the Hazmat/Hospital Sked Stretcher. It features four warm Papoose like pockets that are adjustable for the baby sizes. There is a little hood over each pocket that keeps the baby's head warm and keeps the baby in the pouch if the head end of the stretcher gets lower than the foot end. The babies are carried in line with the Kiddie Litter rather than across to allow for different size babies and to prevent neck injuries during evacuation.

## Kiddie Litter Insert Instruction:



1. The Kiddie Litter insert is attached using straps and buckles on each of its' four corners. The straps are passed through existing slots of the HMH Sked Stretcher and buckled securely as shown in the below illustration.



SAVING LIVES SINCE 1981

VISIT US AT **SKEDCO.COM**  
**1-800-770-SKED**



2. The pouches are opened and the babies wrapped in their blankets are inserted gently into the individual pouches.
3. For smaller babies extra bedding is placed in the pouch to shorten the depth of the pouch.



4. Fasten the velcro strap snugly around the baby being careful not to over tighten it.




5. Pull the little cap over the baby's head.




6. Once all four pockets are utilized fasten the cross straps of the HMH Sked Stretcher and tighten them enough to create plenty of space between the babies and the straps.



7. Pull the foot end of the HMH Sked up and fasten and adjust the foot end straps.



8. There is sufficient room in the foot end of the HMH Sked Stretcher to carry extra diapers and bedding which should be taken with the babies.



9. Be sure to take the baby identification card from the bassinet and insert it into the little window on the corresponding pouch that secures the baby.
10. Once out of the building the stretcher and babies can be covered with a blanket to keep the babies warm and safe.
11. You will notice space between the foam insert and the HMH Sked Stretcher. It is to provide shock absorption when dragging or carrying the HMH Sked/Kiddie Litter thus keeping the babies safer from injuries during the evacuation.





# Skedco® Infant Rapid Access System

The Infant Rapid Access System contains a hanging bag with up to 3 Kiddie Litter Inserts with HMH Skeds. Just pull the looped cord attached to two large zippers and release the Velcro attachment at the top of the bag for access to the Skeds. Perfect for Hospitals evacuation of infants. With 3 Kiddie Litter Inserts with HMH Sked in each Infant RAS, 12 infants can be evacuated at once per system. Custom configurations containing different quantities is available by special order.

## Product Number:

- Infant Rapid Access System with  
(1) Kiddie Litter Insert with HMH Sked: SK-2550-1
- Infant Rapid Access System with  
(2) Kiddie Litter Insert with HMH Sked: SK-2550-2
- Infant Rapid Access System with  
(3) Kiddie Litter Insert with HMH Sked: SK-2550-3



## Infant Rapid Access System Installation



**\*\* Installation of this wall bag should be at least 8 feet (length of the bag) off the floor.**



1. Install screw hooks at joint between wall and ceiling. Leave enough space for grommet to fit onto screw hook. Install one for each grommet in bag inline with grommet holes.

2. Hang wall bag onto screw hooks and store HMH Skeds in bag.



SAVING LIVES SINCE 1981

VISIT US AT **SKEDCO.COM**  
**1-800-770-SKED**